

SHIPMENT OF SAMPLES TO MDS PHARMA SERVICES

- * Ship tube(s) in a leak-proof plastic bag.
- * Add adequate dry ice (20-25 pounds).
- * Complete a sample inventory sheet or write a cover letter with any pertinent data relative to the shipment (i.e. number of samples, type of samples, phase, time points, etc.). This letter or inventory sheet should be placed inside the cooler or between the cooler lid and the shipping box.
- * Pack in a styrofoam cooler that is taped securely, then place into a securely taped fiberboard box.
- * Attach a Class 9 sticker (white and black diamond shape, that can be obtained from your carrier or MDS Pharma Services). Write UN1845 on the outside front of the box next to the Class 9 sticker. Also list the amount of dry ice enclosed in kilograms.
- * Call or fax the Sample Receipt Department at MDS Pharma Services when the shipment has been made and give the appropriate airbill number. Fax (402) 476-7598 - Phone (402) 437-4849.
- * Shipment should be made on Monday, Tuesday, or Wednesday to avoid delivery on weekends. The same care should be taken so the samples are not shipped to arrive on holidays.
- * Ship for next-day delivery to:

MDS Pharma Services
Sample Receipt Department
621 Rose Street
Lincoln, NE 68502
(402) 437-4849

MDS Pharma Services will confirm the receipt of the samples to the sponsor with a follow-up fax.